oxiplexts200[™]

Quantitative Tissue Oximeter









000000

The OxiplexTS200 is a non-invasive, near infrared, dual channel measurement, quantitative tissue oximeter ergonomically designed for clinical use in the OR, ER and doctor's office.





SPECIFICATIONS

Method of Operation	Frequency domain-multiple distance
Modulation Frequency	110 MHz
Measurements	Tissue oxygen saturation
	Oxy- and deoxy-hemoglobin concentration
	Total hemoglobin concentration
	Absorption coefficient
	Reduced scattering coefficient
	Intensity and phase
Light Sources	8 laser diodes emitting at 690 nm
	8 laser diodes emitting at 830 nm
	Laser diodes are time multiplexed
Light Detectors	Photomultiplier tubes
	Computer-controlled bias voltage(Gain)
	Automatic safety shutdown
Average Optical Power	Less than 1 mW
Measurements Channels	Two
Sensors	All fiber optics sensors
	Several sensor types for different applications
	Fiber length up to 10 m
	MRI compatible sensors available upon request
Spatial Resolution	4 Emitter-Detector distances per sensor
Data Acquisition Rate	From 20 ms to minutes
Maximum Experiment Duration	Up to several days (250,000 points)
Software	Windows 8 OS
Electrical Requirements	110-240 Volt, 50/60 Hz
Dimensions	38 cm x 32 cm x 20 cm
Weight	11 kg



www.iss.com

focus and discover

©2016 ISS, INC. ALL RIGHTS RESERVED. All other trademarks or registered trademarks are the property of their respective owners.Caution: Investigational device. Limited by Federal (or United States) law to investigational use. The ISS OxiplexTS200 is presently used for research only. Information furnished by ISS is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or emissions. ISS reserves the right to change the design, specifications, etc., of the product at any time without notice. OxiplexTS200 is covered by US Patents numbers: 5,212,386; 5,492,118; 5,497,769; 5,772,587; 6,078,833; 6,192,261B1. Other US and foreign patents pending.